




# Surgical Guide



**Accurate & precise & High transparency**

Low shrinkage for surgical guide tube tuning. Guides will show high transparency during surgery.



**High mechanical properties**

Strong mechanical properties that resist fracture and distortion.

Designed for the fast, highly accurate production of surgical guides. Possessing industry-leading flexural stiffness and strength, DSG is designed for precise orthodontic work that requires tailored angle and depth.

High flexural strength.

Resists flex and breakage during use.

Nontoxic and biocompatible.



**Strong**



**Transparent**



**Accurate**



**Low Shrinkage**



**DSG**

[ 1kg Bottle ]

**Indications for Use :** The DSG material is indicated for the fabrication of surgical guide, including tissue-supported guides, bone-supported guides, bone reduction guides, and CLP guides.

**Product Description :** DSG is a light-curing resin for the 3D printing of biocompatible dental device using LCD printers (385-405nm).