



Night Guard




Super Clear with LEAP® Technology



High Flexibility & Toughness

DNG enables the direct printing of retainers, nightguards, and dental splints with a clear material and minimal postprocessing. Uniquely engineered for 3D printing at high speeds to enable maximum throughput. DNG possesses industry-leading toughness and biocompatibility. Launching exclusively with our Lux 3 printer.



DNG



1kg Bottle

Applications:

Night guard, Splint, Retainer for aligner, Retainer, Maintainer and other orthodontics appliance.

Indications for Use: The DNG material is indicated for the fabrication of night guards, splints, retainers for aligners, retainers, maintainers and other orthodontics appliance.

Product Description: DNG is a light-curing resin for the 3D printing of biocompatible dental devices using DLP printers (385-405nm).